

MET/TEAM 3.4.0 HOTFIX2 README

This document includes information related to MET/TEAM 3.4.0 HotFix2.

This HotFix resolves various issues in MET/TEAM version 3.4.0. This HotFix is cumulative and contains the changes made in the previous HotFix for MET/TEAM version 3.4.0. See the table below for a complete list.

This HotFix applies only to MET/TEAM version 3.4.0 build 66 or later. Do not attempt to apply this HotFix to any other version of MET/TEAM!

APPLYING HOTFIX TO MET/TEAM AND MET/TEAM CUSTOMER PORTAL

Follow the procedure below to apply this HotFix to MET/TEAM and MET/TEAM Customer Portal. You will need to have Administrator permissions on the machine.

1. Download and extract the contents of the HotFix ZIP archive to a temporary folder on the server where MET/TEAM is deployed.
2. Right-click the ZIP archive and select **Properties** from the popup menu. If there is a checkbox labeled **Unblock** on the lower right of the General tab, check this option and click OK. Otherwise, just click OK.
3. Open the **Internet Information Services (IIS) Manager**
 - a. In the **Connections** pane on the left, expand the **Sites** node and select the **METTEAM** website node.
 - b. In the **Actions** pane on the right, click the **Stop** link in the **Manage Website** section.
 - c. In the **Connections** pane on the left, select the **Application Pools** node.
 - d. In the center pane, select the **METTEAM** item
 - e. In the **Actions** pane on the right, click the **Stop** link in the **Application Pool Tasks** section.
 - f. Repeat these steps for the **METTEAM Customer Portal** application pool.
4. Open **File Explorer**
 - a. Copy the contents of the **\bin** folder that were extracted from the ZIP archive to the MET/TEAM website's **\bin** folder (typically found at C:\inetpub\wwwroot\METTEAM\...), replacing existing file(s).
 - b. Copy the contents of the **\Views** folder that were extracted from the ZIP archive, including all sub-directories, to the MET/TEAM website's **\Views** folder, replacing existing file(s).
 - c. Repeat the above steps for **MET/TEAM Customer Portal** (typically found at C:\inetpub\wwwroot\CustomerPortal\...).
5. Using **IIS**
 - a. In the **Connections** pane on the left, select the **Application Pools** node.
 - b. In the center pane, select the **METTEAM** item
 - c. In the **Actions** pane on the right, click the **Start** link in the **Application Pool Tasks** section.
 - d. In the **Connections** pane on the left, select the **METTEAM** website node.
 - e. In the **Actions** pane on the right, click the **Start** link in the **Manage Website** section.
 - f. In the **Connections** pane on the left, select the **METTEAM** website node.

- g. In the **Actions** pane on the right, click the link in the **Browse Website** section to open MET/TEAM in a browser window
- h. Log in to MET/TEAM
- i. Select the **Help > About** menu and verify the version in the lower left corner of the **About MET/TEAM** dialog indicates version **3.4.0.70**.
- j. Repeat steps 4a to 4h for **MET/TEAM Customer Portal**.

You may now resume using MET/TEAM and MET/TEAM Customer Portal.

Also provided in the ZIP archive is a PDF of the MET/TEAM Help content. Changes have been made to this document as needed that correspond to the tickets indicated in the table below.

MET/TEAM 3.4.0 HOT FIX CHANGES

Key	Issue Type	Description
MTM-15145	FIXED	Resolved an issue that inadvertently prevented results from being added to a Work Order when using MET/TEAM Mobile unless the Work Order was created prior to Mobile Check Out. Also made a change to display the Find Procedure screen when clicking the Add Result toolbar button when the Procedure Used field is blank.
MTM-15205 MTM-15265 MTM-15325 MTM-15408	FIXED	Resolved an issue that may cause unsaved changes to a Work Order to be discarded if an action such as adding, editing, or deleting Parts, Labor, Standards, Accreditations, or Files was performed.
MTM-15235	FIXED	Added Manual Template back in as a Link To option on applicable Find screens.
MTM-15275	FIXED	Resolved an issue that did not populate Model.ToolAssignment.nToolAssignmentUID with the correct value for DataChecks for Tool Assignments during Edit and Check In modes. Also updated documentation for DataCheck.GetExValue() and DataCheck.SetExValue() functions to clarify proper usage and behavior regarding Tool Assignments.
MTM-15295	FIXED	Resolved an issue that caused the nLastStandardCallSheetUID column in the CallSheetStandardLink table to not get populated when adding or swapping a Standard on a Work Order.

Key	Issue Type	Description
MTM-15305	FIXED	Resolved an issue where the Edit Work Order screen would reload when editing a Standard, Accreditation or Part on a Work Order. Also resolved an issue where the Edit Work Order screen would reload when adding results using a Manual Template, causing the Edit Work Order screen to not refresh properly when saving the Manual Template results.
MTM-15335	FIXED	Moved the Bulk Receiving menu item in the Workflow menu down to the second group of items, above Shipping.
MTM-15355	FIXED	Resolved an issue that would cause the server connection to time out when saving a Manual Template with a large number of rows (typically > 100).
MTM-15388	FIXED	Resolved an issue that would cause the prompt indicating signatures will be removed from a Work Order to be displayed twice when performing certain actions on a closed Work Order that has been unlocked for editing.
MTM-15458	FIXED	Minimized reloading the Edit Work Order screen following actions that do not require the Work Order to be saved so that it only happens in situations where digital signatures were removed.