



General Exam Stethoscope

PROFESSIONAL GRADE

PROFESSIONAL DUAL-HEAD STETHOSCOPE

- Better scope. Better care. A quality stethoscope is vital to your success
- Superior Response. Professional Grade. Grace under pressure.

OUTSTANDING ACOUSTIC RESPONSE


- Combination Diaphragm & Bell
- Deliver Wide-Frequency [High/Low] Response
- Proprietary Design Enhances Aural Performance by 50% Compared to Traditional Models

COMFORTABLE TWO-SIDED CHESTPIECE

- Premium Dual-Headed Device Constructed of Durable, Lightweight Aluminum
- Excellent Sound Conduction, Reliable Pressure Response & Quick, Easy Maneuvering
- Non-Chill Bell & Diaphragm Rims Provide Dependable Comfort Against Patients' Skin

FLEXIBLE, LIGHTWEIGHT DESIGN

- High-Quality Headset Hangs Comfortably on the Neck
- Adjustable Aluminum Binaurals
- Angled PVC Ear Tips
- Rotating Stem
- Flexible 32.5" Phthalate-Free Tubing
- Insert Earpiece & Rotate Chestpiece 180° to Change Active Sides



AVAILABLE COLORS

786401
Stethoscope
Dual-Head
Black
50/Case



786402
Stethoscope
Dual-Head
Blue
50/Case



786403
Stethoscope
Dual-Head
Pink
50/Case



Dual-Head
Chestpiece



Wide-Range
Acoustics



Flex-Fit
Headset



Lightweight
Comfort



Versatile
Applications