



AE-711R *Enhancing the Patient Experience*

Headphone

SPEAKER SIZE	27mm
MAGNET MATERIAL	Ferrite
IMPEDANCE	32 ohms @ 500 Hz
SENSITIVITY	103 dB +/-3dB @ 1000Hz
FREQUENCY RESPONSE	20-20000 Hz
MAX INPUT POWER	3mW
CHANNEL BALANCE	Within 6dB
PLUG TYPE	Single 3.5mm TRS
CORD	6' TPE with strain relief
EARPADS	Foam
EAR PIECES	Fixed
HEADBAND	Adjustable, Plastic
WEIGHT	50g
PACKING	100 pieces per carton



ITEM#
1AE6-711RSR-M32

Features



Available in



Perfect Lightweight Personal Headphones

A very affordable entry-level headphone, the 711R is a time-proven design. Compatible with 3.5mm input jack devices, this headphone is the perfect low-cost choice. Providing an adjustable headband, a 6-foot cord that gives a wide range of motion, rotating earcups, and foam ear pads.

Benefits

- Easy to Clean and Store
- Comfortable Experience
- Adjustable Fit Headband
- Supports Rest and Relaxation
- Rotating Ear-cups
- Affordable and Dependable
- Quality Sound
- Improve HCAHP Scores

COMPLIANCE

CPSIA

TPCH

California
Prop 65

