

STERILE GAUZE PADS

WOVEN, 12-PLY, 1's

PRODUCT DESCRIPTION:

- Designed with non-stick material, the gauze pads are comfortable and painless while removing or switching. Made of high grade USP Type VII cotton gauze dressing to withstand both wet-to-moist and wet-to-dry dressing. The soft cotton is ideal for use those with delicate skin types.

FEATURES:

- Highly absorbent, 100% cotton gauze
- Folded edges minimize loose threads and lint
- Packaged in peel-down pouches for convenient, aseptic handling
- Packed in open dispenser boxes for convenient handling



ITEM #	ITEM	QTY
783012	Gauze Pad Woven Sterile 2" x 2" 12-Ply 1's	100/Bx
783013	Gauze Pad Woven Sterile 3" x 3" 12-Ply 1's	100/Bx
783014	Gauze Pad Woven Sterile 4" x 4" 12-Ply 1's	100/Bx

