



STERILE SINGLE USE Transparent Dressing

INDIVIDUALLY WRAPPED DRESSINGS

- Sterile, Single Use Dressings
- Easy two step application
- Keeps contamination out
- Allows for moisture vapor transfer
- Excellent for IV sites, minor burns, and wounds.

FLEXIBLE, LIGHTWEIGHT DESIGN

- Micro porous transparent film with high Moisture Vapor Transmission Rate (MVTR)
- Protects against fluid contaminants & bacteria
- Rounded edges prevent roll up
- Not made with natural rubber latex



**Flexible
Dressing**



**Durable
Construction**



**Maximum
Comfort**



- **ITEM #** 783620
- **DESCRIPTION:**
Dressing
Transparent
Semi-Permeable
- **SIZE:** 4" x 4 3/4"
- **QTY:** 50/Box