



HCG URINE CASSETTE

PREGNANCY TEST

ITEM # 788525

DESCRIPTION:

Pregnancy hCG Test Cassette

QTY: 25/Box

PRODUCT DESCRIPTION:

- The Dealmed hCG Rapid Test Device (Urine) is a rapid visual immunoassay for the qualitative, presumptive detection of human chorionic gonadotropin in human urine specimens. This kit is intended for using as an aid in the early detection of pregnancy.

FEATURES:

- 25 hCG Pregnancy Urine Test Cassettes Individually Wrapped with Disposable Pipettes
- Utilizing a pipette to accurately add the appropriate amount of urine to the test, it can easily confirm a positive test quickly and with complete accuracy.
- Results in 3 Steps: (1) Remove testing device and dropper from pouch. (2) Add 3 drops of urine specimen to the sample well. (3) Read results
- Clear Results: Two Lines = Positive (Pregnant); One Line = Negative (Not Pregnant).



ITEM #	ITEM DESCRIPTION	QTY
788525	Pregnancy hCG Test Cassette	25/Bx



800-569-0570 • sales@dealmed.com • www.dealmed.com