







PADDED ALUMINUM ROLLED SPLINT

ITEM # 787240



DESCRIPTION:

- Rolled Splint, 36"

DIMENSIONS:

- 4.25" W x 36" L

FEATURES:

- Two layers of foam
- Moldable aluminum
- Wipe down exterior

PRODUCT DESCRIPTION:

- Ability to bend into curve shape.
- Moldable foam. Strong aluminum layers.
- Ideal for supportive post injury needs.
- Bend back and remold for a different injury.
- Small, rolled up design allows easy storing.

ITEM #	ITEM	QTY
787240	Rolled Splint, 36"	Ea.

