

Top 5 Things to Know Before Purchasing an ADA-Compliant Wheelchair Scale

Choose a scale that supports patient dignity, staff safety, and compliance confidence.



Space Requirements

Confirm platform clearance and total footprint before you buy. ADA minimum platform space is 36" W x 40" D (32" width may apply if platform height is under 4"). Many scales require a larger footprint for ramps and safe entry/exit.



Third-Party Certification

"ADA-friendly" isn't the same as compliant. Choose independently certified equipment to support accessibility expectations and reduce regulatory risk.



Manufacturer Expertise

Partner with an ADA compliance expert—not a general supplier. The right team helps evaluate space needs, select the right scale, and provide documentation for audits.



Capacity + Accuracy + Tare

Look for clinical-grade accuracy, wheelchair tare, and high capacity (including bariatric patients). Accurate weights support dosing and treatment planning.



Workflow Fit

Choose a scale that supports quick weigh-ins, standardizes processes across staff, and reduces bottlenecks during intake, triage, or rooming.

Additional Clearance

Many scales require extra space for ramps and safe access




Choose Certified ADA-Compliant Scales



Certified compliance protects patients and your organization.

Tight on Space?

Consider a wall-mounted fold-down scale



Top 5 Things to Know Before Purchasing an ADA-Compliant Wheelchair Scale – Why Now Is the Right Time

Health o meter®
Professional


**First to Comply.
Only to Certify.**



First to Align with U.S. Access Board Standards—
Only with Certified ADA Compliant Scales



2920KL-AM



2620KL



2101KL



1100KL

A Note About Health o meter® Professional Scales

As accessibility expectations rise, simply claiming “ADA compliant” is no longer enough. Certified ADA compliance is the difference.

Health o meter® Professional was the first manufacturer to design medical scales aligned with ADA standards and the first to have those designs independently certified.

Today, Health o meter® Professional remains the only manufacturer offering a full portfolio of Certified ADA Compliant scales, supported by documentation facilities can rely on during inspections, audits, and accreditation reviews.

**Don't risk non-compliance.
Choose the only certified
ADA Compliant scales on the
market—trusted by healthcare
providers nationwide.**



Scan to Learn More About Certified ADA Compliant Scales
or visit www.homscales.com/ada-mde