

Accurette™

Mouthpiece for ndd® EasyOne®, Easy on-PC®, EasyOne Pro®, & EasyOne Pro Lab®

The Accurette is the USA-engineered replacement for the ndd® Spirette™

Accuracy confirmed by ndd® spirometer calibration check

Using the SDI Accurette will not affect your spirometer warranty!



Features

- Lower price
- Accurate test results
- Prevents cross-contamination
- No adapters needed

This product is compatible with one of the following: EasyOne®, Easy on-PC®, EasyOne Pro®, and EasyOne Pro LAB®. All brand names are the property of their respective owners. SDI Diagnostics is not affiliated with or endorsed by these brand names or trademarks used solely for compatibility illustration.

Part No.	Description	Qty
29-8060-050	Accurette Mouthpiece	50
29-8060-200	Accurette Mouthpiece	200
29-3150-050	PulmoGuard N with Accurette kit (1 ea)	50

SDI Diagnostics, Inc is your trusted Spirometry Source proudly engineering and precision manufacturing spirometry supplies in the United States for over 40 years

*Accurette*TM

Accuracy Confirmed by ndd[®] Calibration Check

EasyOne[™] DIAGNOSTIC 6.7
© ndd 2000-2010
RecNo 11770

Equipment Information

Software Version 06.07.00.00
Syringe Volume 3.000L

Location _____

Syringe Model _____

Syringe S/N _____

Test Information

Tech ID _____

Calibration Check Results

	Actual	%Diff
Exp. Volume[L]	3.06	1.87
Insp. Volume[L]	2.97	-1.02
Mean Exp. Flow[L/s]	4.13	
Mean Insp. Flow[L/s]	4.37	

Accuracy Confirmed

Equipment Information

Software Version 06.07.00.00
Syringe Volume 3.000L

Location _____

Syringe Model _____

Syringe S/N _____

Calibration Check Results

	Cal1	Cal2	Cal3
Exp. Volume[L]	2.91	2.99	3.04
Mean Exp. Flow[L/s]	0.53	1.16	2.77
%Difference	-2.96	-0.27	1.44

Accuracy Confirmed



This product is compatible with one of the following: EasyOne[®], Easy on-PC[®], EasyOne Pro[®], and EasyOne Pro LAB[®]. All brand names are the property of their respective owners. SDI Diagnostics is not affiliated with or endorsed by these brand names or trademarks used solely for compatibility illustration.

SDI Diagnostics, Inc is your trusted Spirometry Source proudly engineering and precision manufacturing spirometry supplies in the United States for over 40 years