

Heart Rate Display, Alkaline Batteries, Pocket-Sized
with interchangeable probes for vascular and fetal applications

INDUSTRY-LEADING

5

YEAR WARRANTY

ALL PARTS & LABOR
ACCIDENTAL BREAKAGE
NO DIAGNOSTIC FEES
NO HIDDEN FEES
MINIMAL DOWN TIME

MADE IN



REAL PEOPLE REAL SUPPORT

NO ROBOTS.
NO RUNAROUND.
JUST REAL HELP.



The DigiDop II handheld Doppler combines advanced technology with ultimate portability. Made with pride in the USA using high-quality materials, this compact diagnostic tool delivers exceptional Digital Sound Clarity™ technology in a pocket-sized design. Its ergonomic form factor makes it perfect for mobile healthcare professionals.

Like all DigiDop products, it comes with an industry-leading 5-year warranty covering all parts and labor, with zero hidden fees or diagnostic charges.

DigiDop II Features

- Made in the USA with high-quality materials
- Comprehensive 5-year warranty coverage
- Large, bright display
- Ultra-compact, pocket-sized design
- Ergonomic shape for comfortable handling
- Operates on 3 AAA alkaline batteries

DD-770 Ordering Info

Model	Description	Pkg.
DD-770-D2	Heart Rate Display, Alkaline Batteries, 2 MHz fetal probe	1/ea
DD-770-D2W	Heart Rate Display, Alkaline Batteries, 2 MHz waterproof fetal probe	1/ea
DD-770-D3	Heart Rate Display, Alkaline Batteries, 3 MHz fetal probe	1/ea
DD-770-D3W	Heart Rate Display, Alkaline Batteries, 3 MHz waterproof fetal probe	1/ea
DD-770-D5	Heart Rate Display, Alkaline Batteries, 5 MHz vascular probe	1/ea
DD-770-D8	Heart Rate Display, Alkaline Batteries, 8 MHz vascular probe	1/ea
DD-770-D5-SD	Heart Rate Display, Alkaline Batteries, 5 MHz SteadyDop vascular probe	1/ea
DD-770-D8-SD	Heart Rate Display, Alkaline Batteries, 8 MHz SteadyDop vascular probe	1/ea



Interchangeable probes



Seeking a rechargeable alternative of the DD-770?
The DD-770R delivers power on demand.

NEW

The **NEW SteadyDop Probe** simplifies vascular assessments with a stethoscope-like grip and optimized design for easy vessel tracking.



NM0153H