



TSK Microcannulas

TSK STERiJECT™ CSH CANNULA



VALUE PRICED CANNULA

- Regular Handling
- Traditional Product Delivery
- Ultra Thin Wall

BLUNT CLOSED-TIP DESIGN

- Reduced Risk of Trauma
- Easier Handling



ULTRA THIN WALL

- 35% Less Extrusion Force Needed
- Smoother and More Controlled Injection



CODE #	DESCRIPTION	DIMENSIONS
TSK2250SP25	25-PACK STERiJECT CSH Microcannulas	22g x 50mm
TSK2538SP25	25-PACK STERiJECT CSH Microcannulas	25g x 38mm
TSK2550SP25	25-PACK STERiJECT CSH Microcannulas	25g x 50mm
TSK2725SP25	25-PACK STERiJECT CSH Microcannulas	27g x 25mm
TSK2738SP25	25-PACK STERiJECT CSH Microcannulas	27g x 38mm
TSK2750SP25	25-PACK STERiJECT CSH Microcannulas	27g x 50mm
TSK3025SP25	25-PACK STERiJECT CSH Microcannulas	30g x 25mm

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com