



NON-SURGICAL REUSABLE ISOLATION GOWNS

ITEM # 782055

DESCRIPTION:

- Non-Surgical Reusable Isolation Gown, Each

PRODUCT DESCRIPTION:

- Cost conscious and environmentally responsible choice that provides protection from the transfer of microorganisms and bodily fluids
- Used in low- or minimal-risk patient isolation situations and that are not intended for use during surgical procedures, invasive procedures, or when there is a medium or high risk of contamination.

DO NOT USE:

- (I) in a surgical setting;
- (II) where significant exposure to liquid bodily or other hazardous fluids may be expected
- (III) in a clinical setting with high risk of infection.

FEATURES:

- Soft, flexible V-Care barrier fabric tested to meet AAMI / ASNSI PB 70 fluid barrier requirements for Level 2.
- Poly/vinyl fabric is environmentally friendly – PROP65-free, phthalate-free.
- Fabric maintains up to 20 times commercial laundering.
- Overlap back with twill tie neck and back waist closures.
- Generously sized with longer sleeve length for comfortable arm movement.
- Reusable gown may be washed and dried – see washing instructions.

ITEM #	SIZE	QTY
782055	Non-Surgical Reusable Isolation Gown	50/case