

Pulmodyne®

GO-PAP™

Independent PEEP with a Low Flow Connection.

**Just set it and forget it, and concentrate on
caring for your patient.**

The GO-PAP is a disposable emergency CPAP device, that runs off of a low flow connection. It operates at 10LPM, a setting that agrees with all regulators. Like the O2-MAX™ this device delivers approximately 30% FiO2, and provides the option of 5, 7.5 and 10cm of PEEP.

The PEEP settings are independent of flow, so you never have to worry about adjusting your flow to maintain your PEEP levels. This eliminates the need to read from confusing charts or monitor manometers.



*Purchase with or
without Nebulizing
Canister*

The GO-PAP System includes:

- GO-PAP Generator
- Oxygen Supply Line
- BiTracED™ Mask and Headstrap
- OmniClip™

Part No.	Product Description	Qty.
313-4602NEA	GO-PAP with BiTrac ED™ with neb - Adult Medium	10/cs
313-4603NEA	GO-PAP with BiTrac ED™ with neb - Adult Large	10/cs
313-4602EA	GO-PAP with BiTrac ED™ - Adult Medium	10/cs
313-4603EA	GO-PAP with BiTrac ED™ - Adult Large	10/cs

BoundTree
medical

800.533.0523 | www.boundtree.com

*For more information contact your dedicated Account Manager
or call 800.533.0523.*