

Sol-Care®

Safety Winged Blood Collection Needle

-  The Sol-Care® Safety Blood Collection Needle is a safety-engineered device with a sliding safety shield that clicks when activated, telling you the needle is safely locked.



**Sol-Care® Safety
Blood Collection Needle with
Luer Adapter**

Features

- Transparent design that allow visualization of the blood flashback and flow
- Single-handed activation option that provides more flexibility to care for the patient and venipuncture site
- Sliding safety shield pushes over the needle and locks with an audible click, protecting against needlestick injury



Sol-Care® Safety Blood Collection Needle with Pre-attached Holder

- 1 Blood flashback and flow are easy to visualize
- 2 Fixed Needle design with ultra low dead space for reducing medication loss
- 3 Available with or without pre-attached holder
- 4 Sliding safety shield

Sol - Care®

Safety Winged Blood Collection Needle

Sol-Care® Safety Winged Blood Collection Needle with Luer Adapter

Ref. Number	Gauge	Needle Length	Tubing Length	Box Qty	Case Qty
110201030030	21G	3/4"	7"	50	200
110201030031	21G	3/4"	12"	50	200
110201030032	23G	3/4"	7"	50	200
110201030033	23G	3/4"	12"	50	200
110201030034	25G	3/4"	7"	50	200
110201030035	25G	3/4"	12"	50	200

Sol-Care® Safety Winged Blood Collection Needle with Pre-attached Holder

Ref. Number	Gauge	Needle Length	Tubing Length	Box Qty	Case Qty
110201030004	21G	3/4"	7"	25	200
110201030002	21G	3/4"	12"	25	200
110201030008	23G	3/4"	7"	25	200
110201030006	23G	3/4"	12"	25	200
110201030012	25G	3/4"	7"	25	200
110201030010	25G	3/4"	12"	25	200

For more information, please contact:

Sol-Millennium Medical Inc.
311 S. Wacker Dr. Suite 4100
Chicago, IL. 60606



833-637-0857



info@solm.com

NACA-BC-403-1222