

Surgical Gloves

High Quality Premium Sterile Glove

STERILE LATEX SURGICAL GLOVES:

- Curved fingers and rolled cuff
- Excellent comfort and fit
- Micro-rough textured surface
- Superior instrument handling and grip
- Fast and easy to donning
- Washed to minimize skin irritants
- Extra supple with great tactile sensitivity

COMPOSITION:

- Natural Rubber Latex
- Sulfur
- Zinc Oxide
- Organic Accelerators
- Antioxidants
- Synthetic Polymer

PRODUCT SPECS:

- Made with Natural Rubber Latex.
- Powder Free. Chlorinated.
- Micro Textured. Polymer Coated.
- Low Proteins
- Meets ASTM D 3577 & 5151 Compliance Standards
- Single Use Only. Hand Specific.
- Extended Cuff. Beaded Cuff.

RECOMMENDED USES:

- General Surgery
- Orthopedic Procedures
- Outpatient Surgery Centers
- Dental Surgical Procedures
- Veterinary Surgery



ITEM #	SIZE	QTY
787955	5.5	50/Box, 6/Case
787960	6	50/Box, 6/Case
787965	6.5	50/Box, 6/Case
787970	7	50/Box, 6/Case
787975	7.5	50/Box, 6/Case
787980	8	50/Box, 6/Case
787985	8.5	50/Box, 6/Case
787990	9	50/Box, 6/Case

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.