

Pregnancy HCG Test Cassette

Provides early detection of the pregnancy hormone**

Features

- Easy-to-use, one-step process
- Convenient thumb rest for ease in handling
- Easy-to-read results in 5 minutes

Easy-to-read results window

- Early Result Pregnancy Test involves just one easy step
- Hold the absorbent in urine stream at least 5 seconds
- Place test cartridge to rest on a level surface
- Wait 5 minutes for results to show

Early testing

The early result pregnancy test can be used as early as 5 days before you expect your period, which is 6 days earlier than waiting to test until you miss your period.

Over 99% accurate from the day of your expected period*

Accurate at detecting typical pregnancy hormone levels from the day of your expected period.


Note: An individual's hormone levels vary.



ITEM #	QUANTITY
788503	3/Box
788506	6/Box

**Easy-to-read
results window**

**Results in
just 5 minutes**



*At detecting typical pregnancy hormone levels, from the day of your expected period. Note that an individual's hormone levels vary.


**The Early Result Pregnancy Test can be used as early as 5 days before you expect your period. Which is 6 days earlier than waiting to test until you miss your period. In testing, with early pregnancy samples, the Early Result Pregnancy Test gave the following results (percentage of samples from pregnant women given a positive test result): 5 days prior: 68%; 4 days prior: 89%; 3 days prior: 97%; 2 days prior: 98%; 1 day prior: >99%; 0 days prior: >99%.

HCG Pregnancy Test Strip

Dealmed HCG Pregnancy Test Strip
is a rapid, one-step test

Features


- Rapid, one-step process
- Easy-to-follow instructions
- Over 99% accurate at detecting typical pregnancy hormone levels
- Easy-to-read results in 3-5 minutes



ITEM #	QUANTITY
788550	50 Strips/Box

Easy-to-read
results in
3-5 minutes


Control Line
Test Line



POSITIVE



NEGATIVE



INVALID



Dealmed Medical Supplies, LLC

3512 Quentin Rd, Suite 200, Brooklyn, NY 11234 | 800-569-0570 | [dealmed.com](https://www.dealmed.com)