

MetriCide™ OPA Plus Test Strips

Frequently Asked Questions

What is changing with MetriCide OPA Plus Test Strips?

Instructions for use and part number have changed. Please read instructions thoroughly prior to first use of these new test strips. Key product improvements include:

- Better readability with bigger pad size on individual test strips.
- Tamper evident packaging to enhance product safety and reliability providing a clear indication if product has been opened or altered before reaching end user customers.
- Less chance of exposure to humidity and other environmental conditions that can affect performance of the strips by reducing number of strips per package.



	Old Test Strips	New Test Strips
SKU	10-602	10-602A
Product Description	MetriCide OPA Plus Test Strips	MetriCide OPA Plus Test Strips (Replace 10-602)
Quantity per unit	100 strips per bottle, 2 bottles per case	60 strips per bottle, 2 bottles per case
Packaging	Standard packaging	Tamper evident packaging
Shelf Life	15 months	11 months

Why are you making this change?

This change is part of our commitment to continuous improvement and innovation. The new test strips provide better readability, improved packaging and maintain the same high standard of accuracy.

Can I still use the old test strips (10-602)?

Yes, previously distributed MetriCide OPA Plus Test Strips (PN 10-602) remain good until expiration.

What if I still have MetriCide OPA Plus Solution bottles that reference the former test strips part number 10-602 on the label?

Both test strips are acceptable to use with MetriCide OPA Plus Solution. 10-602A will replace 10-602, and the MetriCide OPA Plus Solution label will be updated. Please keep a copy of the Metrex product update memo on file as a reference for your facility.

When will the new test strips be available?

Please contact your local dealer for more details about availability.

Where can I find more information about the new test strips?

For more information, you can download the MetriCide OPA Plus sales sheet, wall charts, and other marketing documents from our website [Metrex.com](https://www.Metrex.com) or contact our customer care team at 800-841-1428 or metrexcustcare@metrex.com.



Scan the QR Code or
visit [metrex.com](https://www.Metrex.com)