




Masimo Sensors RD SET™

Fold-over Style Sensors



> SENSOR CONNECTOR

- > Lightweight
- > No moving parts
- > Smooth edges
- > Sensor insertion keyed
- > Friction fit
- > Tactile and audible feedback at connection
- > Protected from fluid ingress
- > Color coded
- > Finger grips on top of sensor connector

Now
with Improved
SpO₂ Accuracy



> Adt

Adult
>30 kg



> Pdt

Pediatric
10-50 kg




Flexible, lightweight,
flat sensor cable

Soft edges

Color code and
parameter indicator

Lot code

Keyed




SpO₂ accuracy has improved in conditions of motion and no motion, providing clinicians with greater confidence when monitoring oxygenation status.

Product name

Emitter icon

Detector icon




Finger placement guide

Low profile detector

Flex cable between
emitter and detector

Low profile emitter

Fingertip placement guide



Masimo

Masimo Sensors RD SET™

Wrap-around Style Sensors



> SENSOR CONNECTOR

- > Lightweight
- > No moving parts
- > Smooth edges
- > Sensor insertion keyed
- > Friction fit
- > Tactile and audible feedback at connection
- > Protected from fluid ingress
- > Color coded
- > Finger grips on top of sensor connector

Now
with Improved
SpO2 Accuracy




> SOFTOUCH SENSORS

- > Designed to be used whenever skin sensitivity issues are a concern, such as with extremely low birth weight infants
- > Incorporate soft foam and VelAid hook and loop attachment wraps that come in a variety of configurations to address a wide range of clinical uses


> Inf
Infant
3-20 kg



> NeoPt
Neonatal
<1 kg



> Neo
Neonatal
<3 kg or > 40kg



> NeoPt-500
Neonatal
<1 kg



Flexible, lightweight, flat sensor cable

Soft edges

Alignment guide

Color code and
parameter indicator

Flex cable between
emitter and detector

Lot code

Low profile detector

Emitter icon

Detector icon

Product name

Replaceable tape

Keyed

Heat sealed
sensor package

Low profile emitter

SpO2 accuracy has improved in conditions of motion and no motion, providing clinicians with greater confidence when monitoring oxygenation status.



Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

© 2018 Masimo. All rights reserved.



60012/PLM-10451B-1118