

ClearLook® Lighted Ear Curette™



Bionix® ClearLook Lighted Ear Curettes are designed to support precise, controlled cerumen and foreign body removal while helping clinicians deliver a safer, more comfortable patient experience.

Why ClearLook Lighted Ear Curettes Matter

- Illumination and magnification provide clear visibility
- Secure finger placement outside the line of sight and a natural, pencil-like grip support confident technique
- Supports safer, gentler removal that enhances patient comfort
- Single-use design promotes infection prevention and patient safety

ClearLook Lighted Ear Curettes are backed by more than 40 years of Bionix ear, nose and throat innovation.

Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)

International: Contact your distributor or call
+1 419.727.8421

**See clearer. Perform confidently.
Treat safer.**



Trusted performance. Proven results.

ClearLook Lighted Ear Curettes represent the next generation of lighted ear curettes, delivering premium performance that elevates the standard of care in every clinical setting.



**Innovative Solutions
for Everyday Problems**

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)
 800.455.5678 (+1 419.727.8426)
 bionix.com
 bionix@bionix.com (international@bionix.com)

ClearLook® Lighted Ear Curette™



Product Offering:

(Each box includes 50 Curettes + 1 Light Source + 1 Magnification Lens)



#3310 – ClearLook® AngleLoop®
Angled 4mm diameter loop, similar to AngleLoop Red Safe Ear Curette™



#3350 – ClearLook® CeraSpoon®
4mm wide tip, similar to CeraSpoon Yellow Safe Ear Curette



#3320 – ClearLook® MicroLoop®
3mm diameter loop, similar to MicroLoop Green Safe Ear Curette



#3360 – ClearLook® InfantScoop®
2mm wide tip, similar to InfantScoop Blue Safe Ear Curette



#3330 – ClearLook® VersaLoop®
3mm wide tip, similar to VersaLoop Purple Safe Ear Curette



#3345 – ClearLook® Variety Pack
An assortment of all 6 styles of ClearLook Lighted Ear Curettes



#3340 – ClearLook® FlexLoop®
4mm wide tip, similar to FlexLoop White Safe Ear Curette

CPT Reimbursement Information:

Cerumen or Foreign Body CPT Codes

69210 – Removal impacted cerumen (separate procedure), one or both ears.

69200 – Removal foreign body from external auditory canal; without general anesthesia.



**HOW TO
USE VIDEO**

Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)

International: Contact your distributor or call +1 419.727.8421



Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

ClearLook® Lighted Ear Curette™ US Patent #s: 12402974, D1036670, D1090838

©2026 Bionix LLC | MKT-0277 | REV. 03



**Innovative Solutions
for Everyday Problems**

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)
800.455.5678 (+1 419.727.8426)
bionix.com
bionix@bionix.com (international@bionix.com)