

## Vacu-Aide® Suction Machine

HCPCS: E0600,E0600



- Smaller footprint
- Lighter weight
- Quieter
- Longer lasting battery (7325P-D and 7325P-D-EXF)
- 27 LPM (free flow) typical\*\*\*

| PRODUCTS    | Description                              |
|-------------|------------------------------------------|
| 7325D-D     | Vacu-Aide w/o battery, w/internal filtr  |
| 7325D-D-EXF | Vacu-Aide w/o battery, w/external filtr  |
| 7325P-D     | Vacu-Aide w/btty,internl filtr,carry bag |
| 7325P-D-EXF | Vacu-Aide w/btty,externl filtr,carry bag |

| SPECIFICATIONS           | Description                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------|
| Warranty                 | 7325P-D and 7325P-D-EXF: Limited 2 years on unit                                                                  |
| Product Weight           | 7325P-D, 7325P-D-EXF: 3.8 lb. not including power supply 7325D-D, 7325D-D-EXF: 3.4 lb. not including power supply |
| Product Dimensions       | 6.9"H x 7.0"W x 7.1"D                                                                                             |
| Color                    | Blue & Grey                                                                                                       |
| Operating Humidity Range | 15 - 93% non-condensing                                                                                           |
| Size                     | 6.9"H x 7.0"W x 7.1"D                                                                                             |
| Vacuum Range             | 50 - 550 mm Hg +/- 5%***                                                                                          |
| Warranty Limited         | limited 1 year on battery 7325D-D and 7325D-D-EXF: Limited 2 years on unit                                        |

| ACCESSORIES  | Description                                  |
|--------------|----------------------------------------------|
| 22330        | Canister Kit,Drive 1/ea, EA                  |
| 18600-FILTER | Hydrophobic Filter ForSuctMachine,3/bx, BX   |
| 18600-KITN   | UniversalSuctionTubing&FilterKit,1/ea, EA    |
| 22330-12     | Cannister Kit, Drive,12/cs, CS               |
| 610-12B      | 800CC DispoableSuctionCanister&Lid,12/CV, CS |
| 610-48BP     | 800CCDisposableSuctionCanister&Lid,48/CV, CS |
| 7305D-608    | BACTERIA FILTER (NON-STERILE) 12EA/PKG, PKG  |
| 7305D-612    | REPLACEMENT TUBING 6-PK, PKG                 |
| 7305D-614    | ElbowConct,Filter, 3.375" tubing, 1ea/ea, EA |
| 7305D-619    | ELBOW (DFT) 6-PK SS, PKG                     |
| 7305D-632    | Int Fltr Disp Canstr PK48, CS                |
| 7305D-633    | Int Filter Disp Canister Comp Kit Assmd, EA  |

**ACCESSORIES****Description**

|                |                                              |
|----------------|----------------------------------------------|
| 7305D-633-25   | Int Fltr Disp Cnstr Comp Kit Assmd Cs 25, CS |
| 7305D-635      | 800ML DISP CONTR FILTER Cart. Asml,12 P, CS  |
| SUCP TUBING 10 | Suction tubing Blue tipped 10",1/EA, EA      |
| SUCP TUBING 72 | Suction Tubing 72",1/EA, EA                  |