

Document Number: MM1627-XXXXXXX 10280

RE
Originator

Date _____

Quality Systems Process Owner

~~Date~~

Manufacturer Representative

Date

Date	Rev.	CAN #	Reason for Change
12/29/2022	0	N/A	Initial issue

Technical Specification

Title : DynaRide S5 Ultra Lite Wheelchair (K5)

Supplier :

Supplier Item Code: 10280 – DynaRide S5 Ultra Lite Wheelchair (K5), 1pc/cs


Restrictions : No changes may be made to the manufacturing process or any materials without written permission from Dynarex Corporation's QA Department.

1.0 Contract Manufacturer:


2.0 Component Descriptions:

2.1 Materials of Construction (BOM and Drawing)

No.	Spec	Size	Qty	Material	Color
1	Handle	3 adjustment positions	2	Aluminum	Black
2	Backrest	18"	1	PP	Black
3	Armrest Assembly Left	Flip Back	1	Assembly	Black
4	Armrest Assembly Right	Flip Back	2	Assembly	Black
5	Large Wheel	24"	2	Assembly	Black
6	Large Wheel Securing Plate	2 Adjustment Holes	2	Assembly	Black
7	Cross Frame Left	18"	1	Aluminum	Black
8	Cross Frame Right	18"	1	Aluminum	Black
9	Brake Left	N/A	1	Aluminum	Black
10	Brake Right	N/A	1	Aluminum	Black
11	Small Caster	5 adjustment positions	2	Assembly	Black
12	Legrest (Optional)	5 adjustment positions	1	Assembly	Black
13	Leg Frame	N/A	1	Aluminum	Black
14	Right Frame	N/A	1	Aluminum	Black
15	Seat Extension Kit	2" Extension	2	Assembly	Black
16	Seat	16"	1	PP	Black
17	Anti-Tipper	3 adjustment position	1	Steel	Black



2.2 Product Dimensions



2.3 Product Photo




2.4 Others (Product function, performance, weight capacity, electric parameter and so on.)

- Weight capacity is 250 lbs.

3.0 Packaging:


3.1 Bag:



Reorder #	Tolerance ($\pm 20\text{mm}$)		
	Length(mm)	Width(mm)	Put-Up Qty
10280	1100	1360	1 pc/bag



3.1 Master Case:



Reorder #	Tolerance $\pm 10\text{mm}$			Quantity
	Length(mm)	Width(mm)	Height(mm)	
10280	800	320	780	1 pc/cs

4.0 Labeling:

4.1 Graphics:

All color graphics shall be submitted to Dynarex for Approval.

4.2 Lot Numbering:

Lot numbering will be 5-digit purchase order issued by Dynarex.

5.0 Risk Analysis:

Vendor to provide Risk Analysis report for this item.

6.0 Testing & Inspection Criteria:

Testing shall be per contract manufacturers final test procedures. Vendor must provide the Certificate of Compliance (COC) / Analysis (COA) for all the items shipped.

7.0 Warranty:

18 months from date of shipment.

End

UNCONTROLLED COPY