

POWDER FREE NITRILE EXAM GLOVES

SINGLE USE ONLY
MEDICAL EXAM GLOVES

FEATURES:

- Medical Grade. Non-Sterile
- Powder-Free. Latex Free.
- Made in Malaysia
- Made of Synthetic Nitrile Rubber (NBR)
- Beaded Cuff & Textured Fingertips

PRODUCT DETAILS:

- Passes ASTM F1671 test standards for resistance to penetration by blood borne pathogens.
- Complies with 21 CFR parts 170-199 for use in contact with food.
- Manufactured in compliance with ISO 9001, ISO 13485, and US FDA 21 CFR part 820
- Follow Standard and Universal Precautions and are least reactive to chemicals and ensure comfort and tactile sensitivity for delicate patient care.



TESTED
FOR USE WITH
CHEMOTHERAPY
DRUGS



PACKAGING:
Soft-Fill Technology
Packaged flat and
dispense with ease

ITEM #	SIZE	CT BOX/ CASE
787310	Small	100 / 10
787320	Medium	100 / 10
787330	Large	100 / 10
787340	X-Large	90 / 10