

MALE URINAL

For incontinence and medical care.

ITEM # 781148



PRODUCT DETAILS:

- Made of high-quality Polyethylene translucent material for a sturdy and smooth finish
- Easy-grip, notched handle
- Top opening measures 2.5 inches wide
- Designed with wide base to prevent tipping
- Holds up to 1000 cc with 25 cc and 1 ounce increment measurements.

**Attached, Spill
Proof Lid and
Graduated Scale
on three sides of
container**



ITEM #	DESCRIPTION	QTY/BX	QTY/CASE
781148	Urinal Male, Ea.	1	48

FEATURES:

- The light weight container makes it ideal for: hospitals, doctors office, home care facilities, and even at home.
- Great to have on hand for bedridden individuals or for ones on-the-go while road tripping or camping.
- Perfect for the bedridden, wheelchair bound, elderly, and for those who are physically-challenged.