

SINGLE USE STRETCHER SHEET

TWO
SIZES:
40" x 72"
40" x 90"



PRODUCT DESCRIPTION:

- Offer multi-purpose use, such as stretchers, exam tables and other furniture that requires protection.
- The tissue/polyethylene fabric lining provide solid protection from any bodily fluid spills or stains and against dirt and grime.
- With 1-ply tissue/poly material, there is added protection without any compromise to the softness and comfort.



ITEM #	ITEM	QTY
784061	Stretcher Sheet 40" x 72" 1-Ply T/P Blue	50/Cs
784062	Stretcher Sheet 40" x 90" 1-Ply T/P Blue	50/Cs