



POWDER FREE NITRILE EXAM GLOVES

SINGLE USE ONLY
MEDICAL EXAM GLOVES

DESCRIPTION: Powder Free Nitrile Exam Gloves

PRODUCT DETAILS:

- Ideal for all purpose medical use. They provide more protection than vinyl gloves and do not cause any allergies associated with latex gloves.
- These gloves follow Standard and Universal Precautions and are least reactive to chemicals and ensure comfort and tactile sensitivity for delicate patient care.

FEATURES:

- Passes ASTM F1671 test standards for resistance to penetration by blood borne pathogens.
- Fingertip texture for grip on glass wear, small objects, instruments and tools.
- Manufactured in an ISO 9001 certified facility.
- Complies with 21 CFR parts 170-199 for use in contact with food.

**BEADED CUFF • LATEX FREE
TEXTURED FINGERTIP**



PACKAGING:

- Dealmed uses Soft-Fill Technology in their packaging. It is a new standard in packaging and dispensing of Nitrile Gloves. Our gloves dispense with ease because they are packaged flat which allows for more efficient and expedient use.

ITEM #	SIZE	CT BOX/ CASE
787350	X-Small	200 / 10
787360	Small	200 / 10
787370	Medium	200 / 10
787380	Large	200 / 10
787390	X-Large	180 / 10



800-569-0570 • sales@dealmed.com • www.dealmed.com