



STERILE ALCOHOL PREP PADS

INDIVIDUALLY WRAPPED


DESCRIPTION: Alcohol Prep Pads, Sterile, Individually Wrapped

PRODUCT DESCRIPTION:

- Antiseptic, latex-free, gamma sterilized wipes
- Ideal for disinfecting surfaces for injections, medical devices such as stethoscopes and other equipment that contacts patients directly.
- Great additions to first-aid kits, clinics, EMTs, trainers, and at home use medical kits, too.

FEATURES:

- Saturated with 70% isopropyl alcohol.
- Fast, effective, and quick drying technology, prep pads kill more than 99.9% of germs.
- Packaged in an easy-to-open container for quick access.
- Great in size and perfect to place in your pocket or bag.
- Each alcohol prep pad is sealed in its own packet that provides an airtight seal.



ITEM #	ITEM	QTY
781010	Alcohol Prep Pad Medium	200/Bx
781020	Alcohol Prep Pad Medium	300/Bx
781040	Alcohol Prep Pad Medium	400/Bx
781011	Alcohol Prep Pad Large	100/Bx

