

How to check your blood sugar level during pregnancy


Checking your blood sugar regularly can help you to understand what makes your numbers go up or down. Ask your provider how many times you should test daily.



1 Wash your hands.



6 Insert the test strip.



2 Twist the lancet to remove the top.



7 Massage your finger to get more blood to the fingertip.



3 Open the lancing device and insert the lancet.



8 Prick finger and press for drop of blood. The side of the finger is less painful.



4 Adjust the lancing device setting to select the needle depth. Depth will depend on condition of the skin. For softer skin, select 1x. For thicker skin, select 5. This is important for both comfort and effectiveness in drawing blood.



9 Wipe the first drop of blood with tissue. Collect the second drop of blood with the tip of the strip and wait for the results.



5 Turn on the glucometer.



10 Log your results and bring them to ALL your clinical appointments.

Target blood sugar ranges are:

- **Less than 95 mg:** when you wake up or before meals.
- **Below 120 mg:** 2 hours after starting your meal.
- **Above 70 mg:** if on insulin.
- Your target range may be different based on your needs. Talk with your provider to determine what is best for you.

If you have any questions or concerns, please contact us at:

Women's Health Center

Christine Gomez
Comprehensive Perinatal Health Worker
(213) 747-5542, ext. 3454

Adult Department

Nancy Guerrero
Health Educator
(213) 747-5542, ext. 3438

For up-to-date hours of operations, please visit eisnerhealth.org



 eisnerhealth.org

 @EisnerHealth

 @eisner_health