

# Proper storage for x-ray sensors 101

Source: Sodium Dental



# 1

Make sure you loosely coil sensor cable. Cable care is important to help extend life of sensor.



# 2

Properly sterilizing your sensor is very important. Be sure to follow manufacturer instructions for cleaning sensor. Before placing sensor into bin gently.



# 3

Make sure to close lid and store in a safe place. Removing your dental sensor between uses is best way to extend sensor life.



# 4

Make sure you safely store and transport your dental x-ray sensor. This bin can be used to carry sensor between operatories.