UC Irvine MRI Safety Questionnaire

Name: D) Department	epartment:		
Email:		extension:				
impla	ants, de		the device if	you er	creening is to ensure your safety since certain nter a strong magnetic field. Please indicate if el free to ask MRI personnel for clarification.	
→ WARNING If you have one or more of the following, approaching the MRI scanner may cause serious harm or even death.			The	→ CAUTION The presence of any of the following may or may not exclude you from working in the MR environment.		
Pleas	se tell th	e MRI personnel immediately.				
			Yes	No		
Yes	No				Any type of prosthesis (eye, heart valve,	
		Aneurysm clip			limb, penile, etc.)	
		Cardiac pacemaker			Eyelid spring or wire	
		Implanted cardioverter defibrillator (ICD)			Metallic stent, filter, or coil	
		Electronic implant or device			Shunt, vascular access port, or central	
		Magnetically activated implant or device			line.	
		Neurostimulation system			Radiation seeds or implants	
		Spinal cord stimulation system			Swan Ganz or thermo-dilution catheter	
		Bone growth or bone fusion stimulator			Medication patch (nicotine, birth control,	
		Cochlear, otologic or other ear implant			nitroglycerine, etc)	
		Insulin or drug infusion pump			Any metallic fragment or foreign body	
		Pregnancy Prior evo injury			(bullet, shrapnel, etc)	
		Prior eye injury			Wire mesh implant	
					Tissue expander (e.g. breast)	
					Surgical staples, clips, or metallic sutures Wound dressing	
Imn	ortant:				Joint replacement (hip, knee, etc)	
Important: The MRI magnet is always on. Before entering the MR					Bone or joint pin, screw, nail, wire, plate,	
environment you must remove all metallic objects.				Ш	etc.	
These include hearing aids, keys, beepers, mobile					Dentures or partial plates	
phones, eyeglasses, hair pins, barrettes, jewelry, body					Tattoo or permanent makeup	
piercing jewelry, watches, safety pins, paperclips, money					Body piercing jewelry	
clips, credit cards, bank cards, magnetic strip cards,			П		Hearing aid	
coins, pens, pocket knives, nail clippers, and tools.			П	П	Other implant:	
In addition, you must screen the animals that you will be imaging for metal objects. Be aware that all equipment you use to take care of the animal, (i.e. anesthesia machine, scissors, hemostats, etc.) must stay behind the 5 gauss line (yellow and black checkered tape on floor) in MR room.					other implant:	
Sign	ed	(Researcher) Date				
OF	FICE US	SE ONLY				
	ewed by Personn					

Date: _____