FERNO°



The Transcend™ Stair Chair with PowerTraxx® provides maximum safety and power when transporting a patient up or down stairs.

Its electronically-controlled motor, direct drive transmission, and powered tracks push and pull patient loads of up to 550 pounds. The chair and track system carry the weight, increasing patient safety and reducing operator fatigue and risk of injury.



with POWERTraxx technology can transport patients up and down an average of 40 flights of stairs with the touch of a button.



limits track speed during descent or ascent to a safe, manageable level.



provides a comfortable ride on almost any surface with 60% less shock than competitive stair chairs.



offers integrated LED scene lights to illuminate surroundings up to 1 flight of stairs with two brightness options.



Motorized tracks

power the chair up and down stairs.



Extending footrest

provides added support and comfort for patients.



Rear lift handles

with imporved ergonomics for lifting or carrying the chair.

Preventative Maintenance Program



Preventative maintenance program - Parts, labor, travel, two annual preventative maintenance inspections, unscheduled service, and product health report. Replacement parts do not include: batteries,** accesories, other disposable or expendable parts, or damage due to misuse.

**Batteries are included if accompanying Battery Plan is purchased.









- Two battery options
 ENDURACharge 36V runs up to
 40 flights of stairs per charge or
 the Milwaukee® 28V lithium-ion
 battery that is commercially
 available and runs up to 30
 flights per charge. Both are
 rechargeable and easy to
 remove.
- All-metal frame construction provides durability and long life
- Smooth 5-inch front swivel wheels and extra tall 8-inch rear wheels with integrated shock absorbing tires that navigate rough terrain and do not collect debris.
- 5-position lift bar and handles reduce the strain of lifting and provide maximum control for operators of various heights (wireless motor controls work in any handle position)
- Built-in electronic speed control limits track speed during descent to a safe, manageable level

- One inch track-to-ground clearance provides smooth rolling over carpet and rough surfaces, as well as easier descent down stairs
- Extending footrest and lower track angle increase patient sense of security
- Patient seating surface is contoured and has removable injection molded seat panels
- Positive locking hinge mechanism ensures chair is locked in place whether deployed for use or folded for carrying and storage
- The PowerTraxx Conversion
 Kit can retrofit powered tracks
 to the manual Transcend Stair
 Chair.
- Easy to operate with pushbutton wireless controller



SPECIFICATIONS

Item Number	Base model: 0731403
Height (min)	36 in / 908 mm
Height (max)	60 in / 1518 mm
Width (overall)	20 in / 505 mm
Width (seat)	16 in / 408 mm
Depth (front to back)	30 in / 760 mm
Length (folded)	36 in / 912 mm
Width (folded)	20 in / 507 mm
Depth (folded)	7.9 in / 200 mm
Weight	51 lb / 23.3 kg
Load Limit	550 lb / 250 kg
Depth (front to back)	30 in / 760 mm
Power Supply	36V or 28V Lithium Ion Battery

TRANSCEND OPTIONS



Wheel Locks



Extending Lift Bar



Footrest



Armrests



Rear Lift Handles



Glide Bar

Headrest



No. 0822023



Rear Lift Handles No. 0822028



Restraint Set No. 0822024



Linkage/Close Handle Kit No. 0822536



PATIENT HANDLING ACCESSORIES

ENDURACharge Charger No. 0822509

ENDURACharge 36V Bat. No. 0822484



Milwaukee Charger No. 0832199 Milwaukee 28V Battery



O2 Holder Kit No. 0822053

No. 0832198

ACCESSORIES NOT PICTURED: PRO F1 Compatibility Kit & Knee Gatch