


**X Series<sup>®</sup> Advanced**

**ZOLL<sup>®</sup>**



**ADVANCE YOUR CARE**

# THE GROUNDBREAKING NEW MONITOR/DEFIBRILLATOR FROM ZOLL



To improve patient outcomes, healthcare providers rely on feedback and data at every step of the rescue. Access to real-time clinical feedback on the scene can take your care to the next level and enable you to manage multiple disease states.

## Real-Time Clinical Feedback—Available Only Through ZOLL

For years, ZOLL® devices have provided real-time clinical feedback on CPR to help rescuers deliver high-quality compressions to patients in cardiac arrest. The X Series® Advanced monitor/defibrillator goes beyond cardiac arrest, providing real-time clinical feedback for patients in respiratory distress and patients with a traumatic brain injury:



**Real CPR Help®** – provides real-time clinical feedback on compression rate and depth



**Real BVM Help™** – provides real-time clinical feedback on ventilation volume and rate



**TBI Dashboard™** – provides the right information at the right time to help rescuers effectively manage traumatic brain injury (TBI) patients

# DISEASE STATE MANAGEMENT AND ACCESS TO KEY DATA




## Manage Multiple Disease States

In the pre-hospital environment, access to data helps you manage trauma, cardiac arrest, chest pain, respiratory distress, pediatric issues, and other conditions. X Series Advanced provides information and feedback on CPR, ventilation, and TBI, increasing the range of treatable patients when compared to other currently available monitor/defibrillators.

## Seamless Device Integration


For cardiac arrest patients, the X Series Advanced features Real CPR Help and Real BVM Help, providing simultaneous real-time feedback on both CPR and ventilation quality. When used with the AutoPulse®, ZOLL's automated CPR device, the X Series Advanced captures the unique chest compression data and displays it on screen.

Standard Monitor/defibrillator



- Cardiac Arrest
- STEMI
- Other

X Series Advanced



- Cardiac Arrest
- STEMI
- Respiratory
- Trauma
- Other

## Easy Access to Data

With RescueNet® CaseReview, clinicians can view data from the X Series Advanced post-case for effective QA/QI—no more chasing crews to get case files. With the push of a button, medics can send files directly to your secure system for performance-focused debriefing. Access to data, including ventilation and CPR performance, is as fast as opening a web browser.





## The Power of Telehealth

With ZOLL's remote viewing telehealth technology, healthcare professionals can view data from the X Series Advanced from a remote location. Remote viewing capability provides effective decision support that enables EMS providers to meet their telehealth

### X Series® Advanced transmission



Cellular tower



Remote view server powered by AWS



Web browser



NEED MORE  
INFORMATION?

Please contact your local ZOLL representative or distributor/partner, or email [ServiceContractsAdmin@zoll.com](mailto:ServiceContractsAdmin@zoll.com)



## ZOLL MEDICAL CORPORATION

An Asahi Kasei Group Company | 269 Mill Road | Chelmsford, MA 01824 | [www.zoll.com](http://www.zoll.com)

Copyright © 2020 ZOLL Medical Corporation. All rights reserved. X Series, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners. ZOLL X Series Advanced is not available for sale in the United States. The product has not received regulatory clearance/approval by the Food and Drug Administration.  
Printed in U.S.A. MCN EP 2010 0329

# ZOLL®